

Odil: implementation of a DICOM library in C++

Julien Lamy

ICube, Université de Strasbourg-CNRS



C++ DICOM libraries today

- We're developing DICOM tools:
 - In C++
 - As free software
 - Using network services, client *and* server
- DICOM libraries (C++, free software):
 - Conquest, DICOM SDL: no API documentation, dead (?)
 - dicom3tools, Imebra: no network support
 - GDCM: incomplete network support (client only)
 - Winner by default: DCMTK

Why develop a new library?

- API is getting old

- Duplication of many standard library features (OFList, OFVector, OFString, OFStream, etc.)
- Errors not based on exception: expressed through OFCondition

```
OFString value;  
OFCondition const condition =  
    dataset.findAndGetOFString(DCM_PatientID, value);
```

- No const-correctness, duplication of code due to lack of templates, callbacks are **void** (*) (**void***), etc.

- Inconsistent memory management

- Some functions take pointers, other take references
- Transfer of ownership not always specified (memory leak)

- Documentation not up to date: need to look in .cpp files, undocumented parameters

A newcomer: Odil

- Started as a set of wrappers around DCMTK: DCMTK++
- Quickly evolved to native implementation
- Time-line:
 - April 2015 (v0.1): DCMTK wrappers (data set accessors, network, association, messages)
 - June 2015 (v0.2): native data set (with conversion to and from DCMTK), JSON
 - September 2015 (v0.3): native reading and writing, native Basic Directory (i.e. DICOMDIR) creation, XML
 - January 2016 (v0.4): rename to Odil, native association, DUL, services (C-STORE, C-FIND, C-MOVE, C-GET)

Language: C++11 goodness

- C++ 11 brings cleaner-looking code, easier to maintain
 - Type inference:
`auto const & x = std::min(v.begin(), v.end());`
 - Range-for: `for(auto const & item: sequence) { /* ... */ }`
 - Lambda expressions:
`auto const f = [] (int x) { return x+1; }`
 - _INITIALIZER list, threading, algorithms, variadic templates, etc.
- Get help from Boost for missing features (XML, networking, regular expression, etc.)
- Bleeding edge not necessary: Odil can be built on Debian 7, Ubuntu 12.04, CentOS 6

Dictionary

- List of known elements: tag, keyword, VR, VM
- Registry: namespace holding the tags
- Dictionary: mapping a tag to its keyword, VR and VM
- Automated generation using the Docbook version of the DICOM standard

```
// Get tag from registry  
auto const patient_name = odil::registry::PatientName;  
// Get information from dictionary  
auto const entry = odil::public_dictionary.find(patient_name).second;  
std::cout  
    << entry.name << ": " << entry.vr << ", " << entry.vm  
    << "\n";
```

Data set

- VR and types: reduce the 31 VRs to basic types (integers, reals, strings, binary)
- Every element is an array of values

```
odil::DataSet data_set;  
data_set.add(odil::registry::PatientName, { "Doe^John" });  
// Tag-based access  
data_set.as_string(odil::registry::PatientName) = { "Doe^Jane" };  
// Keyword-based access (slower)  
std::cout << data_set.as_string("PatientName", 0) << "\n";
```

I/O: binary

- Stream-based (std or boost::filesystem)
- Field-tested with with various transfer syntaxes (little endian with implicit and explicit VR, big endian, JPEG lossless, etc.)
- Choice of transfer syntax, partial reading, optional group length, etc.

```
auto const header_and_data_set = odil::Reader::read_file(  
    std::ifstream("foo.dcm"));  
std::cout <<  
    header_and_data_set.second.as_string("PatientName", 0) << "\n";  
odil::DataSet header;  
header.add("SourceApplicationEntityTitle", {"MYSELF"});  
odil::Writer::write_file(  
    header_and_data_set.second, header, std::ofstream("bar.dcm"));
```


I/O: JSON, XML

- Standard representation of data sets in JSON and XML
- Both input and output
- JSON: JSONCpp, XML: Boost.ptree

```
std::ifstream istream("foo.json");  
Json::Value json;  
istream >> json;
```

```
auto const data_set = odil::as_dataset(json);  
std::cout << data_set.second.as_string("PatientName", 0) << "\n";
```

```
auto const xml = odil::as_xml(data_set);  
std::ostream ostream("foo.xml");  
boost::property_tree::xml_parser::write_xml(ostream, xml);
```

Networking: association

- Under the hood: Boost.Asio
- Handles both TCP/IP and DICOM layers

```
odil::Association association;  
association.set_peer_host("184.73.255.26");  
association.set_peer_port(11112);  
association.update_parameters()  
    .set_calling_ae_title("myself")  
    .set_called_ae_title("AWSPIXELMEDPUB")  
    .set_presentation_contexts({  
        {  
            1, odil::registry::PatientRootQueryRetrieveInformationModelFIND,  
            { odil::registry::ExplicitVRLittleEndian }, true, false  
        }  
    })  
});
```

Networking: services

- “Classic” services (C-STORE, C-FIND, C-GET, C-MOVE)
- Callback-based (function, functor or lambda)

```
odil::DataSet query;  
query.add("PatientName", { "Doe^*" });  
query.add("PatientBirthDate");
```

```
odil::FindSCU scu(association);  
scu.set_affected_sop_class(  
    odil::registry::PatientRootQueryRetrieveInformationModelFIND);  
scu.find(  
    query,  
    [](odil::DataSet const & data_set)  
    {  
        std::cout << data_set.as_string("PatientName", 0) << ": "  
            << data_set.as_string("PatientBirthDate", 0);  
    });
```

Networking: services

■ SCP example

```
odil::Association association;  
odil::EchoSCP echo_scp(association,  
    [](odil::message::CEchoRequest const &)  
    {  
        std::cout << "Somebody pinged us!\n";  
        return odil::message::Response::Success;  
    });  
  
association.receive_association(boost::asio::ip::tcp::v4(), 11112);  
auto const message = association.receive_message();  
echo_scp(message);
```

Availability

- Free software, examples available with the source code: everybody welcome aboard!
- Works on Linux and OS X
- GitHub repository: <https://github.com/lamyj/odil>
- Packaged officially in Debian (testing, unstable) and Ubuntu (xenial), thanks to Debian-Med
- Unofficial packages: <https://github.com/lamyj/packages>
- Used in other projects
 - Dicomifier: DICOM converter
 - Dopamine: document-oriented PACS

Work in progress

- Python wrappers

- Data set and I/O: done

- ```
header, data_set = odil.read("foo.dcm")
print data_set.as_string("PatientName")[0]
```

- Networking, XML, JSON: in progress

- Web services (STOW, QIDO, WADO)

- OS X packaging (Homebrew), Windows